

DR Upgrade – choice without compromise.



Konica Minolta quality and reliability without a connection to the generator.

Konica Minolta has created PerformanceDR, the wireless DR solution that offers the convenience of Automatic Exposure Detection (AeroSync®) eliminating the connection to the generator, and the outstanding performance of Konica Minolta Digital Radiography detectors.

- AeroSync - No connection to the generator needed
- Choice of ULTRA or CS-7 image acquisition
- Industry leading 7-year warranty on the power cell



Just what you need.

Superior image quality

CS-7 and Ultra feature Konica Minolta's Hybrid (H) Processing; a unique image processing technique developed to realize superior image quality and improved patient diagnosis and easily integrate with existing RIS and PACS infrastructure to enable the technologist to acquire, view, and complete a study in a smooth, time saving workflow.

The right size for your practice

Private Practice and Outpatient Clinics requiring wireless DR productivity and that have already implemented PACS. It includes Konica Minolta's PerformanceDR flat panel detector and Ultra Software that simplifies DR image acquisition. PerformanceDR offers only the features needed for primary acquisition.

Key Features		
PerformanceDR	ULTRA OPERATING PLATFORM	CS-7 OPERATING PLATFORM
<ul style="list-style-type: none">• Up to 125 exposures/ 3.4 hours• ~ 10 second Panel Refresh Time (SSRM ~ 6 seconds) (Cycle Time)• Meets International Specifications for true 14x17 cassette size• Wireless Connectivity Package (includes Access Point)• Lithium Ion capacitor includes exclusive 7 year Warranty against replacement for reduced operational cost• AeroDR Docking Station<ul style="list-style-type: none">• FPD Activation Module• FPD Recovery Module• AeroSync Connectivity	<ul style="list-style-type: none">• Study Append• Free Text Annotation• Automatic Masking• Study List Filter• History Search• Image Zoom• Grid Suppression• HIPAA compliance enabling features (audit trail, auto log-out)• Procedure Code Mapping Tool• Foldering• DICOM Store• DICOM Modality Worklist• Image Gently / Image Wisely• Image Merge	<ul style="list-style-type: none">• Study Append• Free Text Annotation• Automatic Masking• Study List Filter• History Search• Image Zoom• Grid Suppression• HIPAA compliance enabling features (password protection, screen black out)• Procedure Code Mapping Tool• Study Combine and Move• DICOM Store• DICOM Modality Worklist

Comprehensive Coverage for Peace of Mind

Blue Moon TCO for PerformanceDR is designed to maximize customer satisfaction throughout the lifecycle of the imaging system while helping customers manage their cost of ownership.

The support approach from Konica Minolta, in conjunction with Authorized Reseller Partners, creates a simple path for ongoing support. Blue Moon support plans help deliver value and confidence while ensuring the worry-free operation of every PerformanceDR system.

Purchase up front as a 5-year TCO option for the lowest long term cost of ownership or pay annually to suit the financing requirement.



Konica Minolta Healthcare America, Inc.
411 Newark Pompton Turnpike
Wayne, New Jersey 07470
Tel: (973) 633-1500 Fax: (973) 523-7408
konicaminolta.com/medicalusa

M1235 1017 RevA